

## De-multiplexing Cell3™Target Data – Illumina Sequencing

This document outlines an example procedure for sample de-multiplexing of pooled Cell3™Target libraries following Illumina sequencing. The i7 indexing read contains a 9bp molecular tag in addition to the unique 8bp sample index. The i5 index contains 8bp unique sample index only. The final read structure is Y\*, I8Y9,I8,Y\*.

### Input requirements

1. Illumina bcl2fastq software – a Base Call (BCL) Files to Fastq conversion software  
[http://emea.support.illumina.com/sequencing/sequencing\\_software/bcl2fastq-conversion-software.html](http://emea.support.illumina.com/sequencing/sequencing_software/bcl2fastq-conversion-software.html)
2. Basecalls directory
3. Samplesheet.csv – ensure there are no N's present in the sample sheet (examples provided on request)
4. Check RunInfo.xml to ensure the correct number of cycles have been performed for the indexes
  - a. i7 = 17 cycles
  - b. i5 = 8 cycles

### Expected file outputs per sample

1. Read 1 fastq file (R1)
2. Read 2 fastq file (R3)
3. fastq file containing information on a per fragment molecular tag (R2)
4. Index 1 and 2 fastq files (these are not required)

### Procedural steps

1. Open a command terminal
2. Move into the basecalls directory of the data to be de-multiplexed
3. Run the following bcl2fastq command  

```
bcl2fastq --create-fastq-for-index-reads --mask-short-adapter-reads 0 --use-bases-mask Y*,I8Y9,I8,Y* --no-lane-splitting (optional)
```

\* replace with the number of cycles performed for read 1 and 2

