

Protocol Guide v1.2.2

Cell3™ Target

DNA Target Enrichment for
Next Generation Sequencing
(Illumina Sequencers)





Table of Contents

Workflow Overview	3
Chapter 1: Library Preparation	4
1.A Kit version (a): End-repair / A-tailing (non-fragmentation)	5
1.B Kit version (b): Fragmentation and end-repair / A-tailing	6
1.C Ligation of Illumina UMI adapters	7-8
1.D Library amplification	9-10
1.E Library QC	10
Chapter 2: Probe Hybridization and Capture Enrichment	11
2.B Probe capture on Streptavidin beads and washes	12-13
2.C Captured library amplification and clean-up	14-15
2.D Captured library QC	15



Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

Workflow Overview

a) For use with cell free DNA as input

b) For use with genomic DNA as input





Figure 1. Flow chart outlining the main steps of the Cell3™ Target workflow. Blue boxes refer to library preparation steps (3h); while green boxes refer to probe hybridization / capture and target enrichment steps (8h).



Chapter 1: Library Preparation

Input DNA requirements


- Determine DNA concentration using a fluorometric method (e.g. Qubit assay, Invitrogen).
- Resuspend in molecular biology grade water or a low EDTA buffer.
- **1-100 ng** DNA is recommended when using UMIs
- **100-1000 ng** DNA is recommended if UMIs are not required or for PCR-free libraries

Input DNA requirements for FFPE samples

The following table provides a guideline on **FFPE DNA** input quantities to use according to the DIN score observed:

Input DNA guidelines for DNA samples extracted from FFPE.			
DNA input parameters	DIN score >8	DIN score 3-8	DIN score <3
Fold increase compared to high-purity DNA	No increase required	Increase input DNA quantity by 1.5-4-fold	Increase input DNA quantity by 5-10-fold

- A minimum input of **10 ng FFPE DNA** is recommended irrespective of the DIN score.



Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

1.A Kit version (a): End-repair / A-tailing (non-fragmentation)

Before you start

- Thaw the End-repair / A-tailing Buffer (10x) (red cap) and Ligation Buffer (5x) (blue cap) at RT and mix by briefly **vortexing**, brief spin, **keep on ice**.
- Mix the End-repair / A-tailing Enzyme Mix (5x) (red cap) and DNA Ligase Enzyme (blue cap) by **tapping the tube**, brief spin, **keep on ice**.

Procedure

1. Set up the ER/AT program.


Step	Temperature	Time
1	4°C	Hold
2	20°C	30 min
3	65°C	30 min
4	4°C	Hold

Set the heated lid at **70°C**, volume to **50 µl** and **start**.

Prepare the reaction, **keep on ice**. Mix by **vortexing or pipetting**, spin.

Components	Volume for 1 reaction
End-repair / A-tailing Buffer (10x)	5 µl
DNA sample	X µl
Nuclease-free water	(35 - X) µl
Total	40 µl

2. Add **10 µl of End-repair / A-tailing Enzyme Mix (5x)** for a total volume of **50 µl**. Mix by **vortexing or pipetting 10-15 times**, spin.
3. Transfer to the pre-chilled thermocycler (4°C) and “skip” to the next step.
4. **Keep on ice. Immediately proceed to the ligation step (1.C).**



Protocol Guide v1.2.2

Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

1.B Kit version (b): Fragmentation and end-repair / A-tailing

- DNA input amounts lower than **50 ng**, add the **Fragmentation Enhancer**.
- For longer insert sizes see **Appendix I** of full protocol.

Before you start

- Thaw the Fragmentation Buffer (10x) (**red** cap), the Ligation Buffer (5x) (**blue** cap) and the Fragmentation Enhancer (**orange** cap) (if required) at **RT** and mix by **briefly vortexing**, spin, keep on ice
- Mix the Fragmentation Enzyme Mix (5x) (**red** cap) and the DNA Ligase Enzyme (**blue** cap) by **tapping** the tube, spin, keep on ice.

Procedure

1. Set up the program in a thermocycler.


Step	Temperature	Time
1	4°C	Hold
2	20°C	30 min
3	65°C	30 min
4	4°C	Hold

Set the heated lid at **70°C**, volume to **50 µl** and **start**.

2. Prepare the reaction, **keep reaction on ice**. Mix well by **vortexing or pipetting**, spin.

Components	Volume for 1 reaction (DNA input ≥ 50 ng)	Volume for 1 reaction (DNA input < 50 ng)
Fragmentation Buffer (10x)	5 µl	5 µl
DNA sample	X µl	X µl
Fragmentation Enhancer	-	2.5 µl
Nuclease-free water	(35 - X) µl	(32.5 - X) µl
Total	40 µl	40 µl

3. Add **10 µl** Fragmentation Enzyme Mix (5x) for a total volume of **50 µl**. Mix by **vortexing or pipetting** **10-15** times, **brief spin**.
4. Transfer to the thermocycler (4°C) "skip" to the next step.
5. **Keep on ice. Immediately proceed to the ligation step (1.C).**



Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

1.C Ligation of Illumina UMI adapters

Before you start

Equilibrate the **Target Pure™ NGS clean-up beads** to room temperature for 20-30 minutes. Remove the required number of Illumina UMI Adapters, **thaw on ice** and spin.

Procedure

1. Set up the program in a thermocycler.


Step	Temperature	Time
1	4°C	Hold
2	20°C	15 min

Set the lid to **“not heated”** (or leave the lid open), volume to **100 µl** and **start**.

2. Illumina UMI Adapters are provided at a concentration of 15 µM.
 - Use at **1.5 µM for DNA inputs less than 50ng** (1:10 dilution)
 - Use directly from the tube for **≥50ng**.
3. While keeping the end-repaired / A-tailed DNA samples **on ice**, add 5 µl of Illumina UMI adapter to each sample and **mix** by **pipetting or briefly vortexing**.
4. Prepare the ligation buffer master mix. Mix by vortexing, spin and **keep on ice**.

Components	Volume for 1 reaction
Ligation Buffer (5x)	20 µl
DNA Ligase Enzyme	10 µl
Nuclease-free Water	15 µl
Total	45 µl

5. Add **45 µl** of ligation master mix for a total volume **100 µl**.
6. Mix by pipetting **10 times or brief vortexing** and spin.
7. Transfer to the pre-chilled thermocycler (4°C) and “skip” to the next step.
8. Remove samples and proceed to the clean-up step using **Target Pure™ NGS clean-up beads**.




Protocol Guide v1.2.2

Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

Clean-up of adapter ligated library

9. Add **90 µl** of **Target Pure™ NGS clean-up beads** to a fresh tube for each sample.
10. Transfer **100 µl** of **adapter ligation reaction** to the 90 µl of Target Pure™ NGS clean-up beads and mix well by pipetting **15-20 times**.
11. Incubate for **5 minutes** at **RT**.
12. Prepare a **fresh** solution of 80% ethanol / 20% molecular biology grade water (800 µl per sample).
13. Place tubes on the magnet for **5 minutes**.
14. Remove and discard the **supernatant**, without disturbing the pelleted beads.
15. Add **200 µl** of 80% ethanol to the tube/well and incubate for **30 seconds**.
16. Repeat steps 14-15 for a total of **two 80% ethanol washes**.
17. Keeping the tubes on the magnet, slowly remove and discard the **supernatant**.
18. Remove any residual liquid from the tubes/wells with a **10 µl** pipette.
19. Keeping the tubes on the magnet, incubate at **RT** with **open lids** for **5 minutes** or **until the beads are dry. Do not over dry**.
20. Remove the tubes from the magnet and resuspend in **27 µl** of **Buffer EB** or **molecular biology grade water**.
21. Incubate the tubes for **2 minutes** at **RT**.
22. Place the tubes on the magnet for **2 minutes** at **RT** to pellet beads.
23. Recover **24 µl** of **supernatant** and transfer it to a fresh **1.5 ml low-bind tube**.

Stop point: adapter ligated libraries can be stored at **-20°C**



Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

1.D Library amplification

Before you start

- Thaw the PCR Master Mix – PreCap (2x) (**green** cap) and the Primer Mix – PreCap (10 µM) (**black** cap) on ice. Mix by **vortexing** and spin. Keep **on ice**.
- Equilibrate the **Target Pure™ NGS clean-up** beads to **RT** for **20-30 minutes**.

Procedure

1. Set up the following program in a thermocycler.


Step	Temperature	Time	Cycles
1	98°C	Hold	1
2	98°C	2 min	1
3	98°C	20 sec	4-12*
4	60°C	30 sec	
5	72°C	30 sec	
6	72°C	1 min	1
7	4°C	Hold	1

* as a general guideline: **4 cycles for 200 ng, 5 cycles for 100 ng, 6 cycles for 50 ng, 9 cycles for 10 ng, 12 cycles for 1 ng.**
 For FFPE DNA increase the cycle number by 1 (e.g. **7 cycles for 50 ng**).

Set the heated lid at **105°C**, volume at **50 µl** and **start**.

Components	Volume for 1 reaction
PCR Master Mix – PreCap (2x)	25 µl
Primer Mix – PreCap (10 µM)	2.5 µl
Total	27.5 µl

2. Prepare the PCR reaction master mix in a separate PCR tube for each sample **on ice**.
3. Mix by pipetting **10 times** or briefly **vortexing** and spin.
4. Transfer **22.5 µl** of adapter-ligated library to the **27.5 µl** of PCR reaction mix. Mix by pipetting **10 times** or brief **vortexing** and spin.
5. Transfer to the pre-heated thermocycler (**98°C**) and **skip** to the next step.
6. Remove the samples from the cycling block and **proceed immediately to library clean-up using Target Pure™ NGS clean-up beads**.



Protocol Guide v1.2.2

Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing



Clean-up of amplified library

7. Add **50 µl** of **Target Pure™ NGS clean-up beads** to a fresh tube for each sample.
8. Transfer the **50 µl** of PCR amplified library to the 50 µl of Target Pure™ NGS clean-up beads and mix well by pipetting **15-20 times**.
9. Incubate for **5 minutes** at RT.
10. Place the tubes on the magnet for **5 minutes** at **RT**.
11. Keeping the tubes on the magnet, remove and **discard** the **supernatant**.
12. Add **200 µl of 80% ethanol** to the tube/well and incubate for **30 seconds**.
13. Repeat steps 10-11 for a total of **two 80% ethanol** washes.
14. Keeping the tubes on the magnet, remove and discard the **supernatant**.
15. Remove any residual liquid from the tubes/wells with a 10 µl pipette.
16. Keeping the tubes on the magnet, incubate at **room temperature, lids open** for **5 minutes** or until the **beads are dry. Do not over dry**.
17. Remove the tubes from the magnet and resuspend in **32.5 µl of nuclease-free water**.
18. Incubate the tubes for **2 minutes** at **RT**.
19. Place the tubes on the **magnet** for **2 minutes** at **RT**.
20. Carefully recover **30 µl of supernatant** and transfer it to a fresh **1.5 ml low-bind tube**.

Stop point: 4°C overnight or at -20°C for long term storage.

1.E Library QC

- DNA concentration (ng/µl) and total yield (ng) using fluorometric method (e.g. dsDNA BR assay kit, Invitrogen **Qubit**)
- DNA quality in terms of expected fragment size distribution and absence of additional lower or higher molecular weight peaks (recommended) **TapeStation**.



Chapter 2: Probe Hybridization and Capture Enrichment

Before you start

1. Switch on a vacuum concentrator and set the temperature to **70°C or lower**.
2. Thaw the Hybridization Buffer (2x) (**blue** cap), the Hybridization Enhancer (**brown** cap), the Universal Blockers (**orange** cap) and the COT-1 Human DNA (**red** cap) at **RT**. Thaw the Cell3™ Target: Probe Set **on ice**.
3. Mix each component by **vortexing** and **spin**.

NOTE: Heat the Hybridization Buffer (2x) (**blue** cap) at **65°C in a heat block** and vortex every few minutes if crystallisation is seen

Procedure

1. Set up the following program on a thermocycler.

Step	Temperature	Time	Cycles
1	95°C	Hold	1
2	95°C	30 sec	1
3	65°C	4 hours	1
4	65°C	Hold	1

Set the lid to heated at **100°C**, the sample volume to **17 µl** and **start**.


2. Pool equal amounts individual libraries into **1.5 ml low-bind** tube to a total of **1000 ng**.
3. Add **5 µl of COT-1 Human DNA** and **2 µl of Universal Blockers** to the library pool. Mix by vortexing and spin.
4. Place the tube (lid open) in the vacuum concentrator until pool is **completely dry**.

Stop point: the dried down library pool can be stored **overnight at 4°C**.

5. Prepare the hybridization reaction mix by adding to the dried down pool.

Components	Volume for 1 reaction
Hybridization Buffer (2x)	8.5 µl
Hybridization Enhancer	2.7 µl
Cell3™ Target: Probe Set	4 µl
Nuclease-free water	1.8 µl
Total	17 µl

6. **Mix** by pipetting **10 times**, briefly centrifuge and incubate at **RT for 10 minutes**.
7. Transfer hybridization reaction mix to a **0.2 ml PCR tube**, brief spin.
8. Place the hybridization reaction mix in the thermocycler (95°C), skip to the next step.
9. Incubate the hybridization reaction mix at **65°C** on the thermocycler for **4-16 hours**.



Protocol Guide v1.2.2

Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

2.B Probe capture on Streptavidin beads and washes

Before you start

- Equilibrate the **Dynabeads® M-270 Streptavidin** to **RT for 30 minutes**.
- Thaw all Wash Buffers (S, 1, 2, 3, B) at **RT, vortex** and spin.
- Heat Wash Buffer 1 at **65°C** to resuspend precipitated particles.

IMPORTANT: for **Cell3™ Target Exome Panel** or **custom probe sets of Tier-4** and above, refer to **Appendix V** to proceed with the alternative post-hybridization capture protocol.

Preparation of wash buffers


1. Dilute for each capture reaction. Prepare a 1x working solution in **1.5 ml** tubes:

Components	Stock solution	Nuclease-free water	Total
Stringent Wash Buffer (10x)	40 µl	360 µl	400 µl
Wash Buffer 1 (10x)	30 µl	270 µl	300 µl
Wash Buffer 2 (10x)	20 µl	180 µl	200 µl
Wash Buffer 3 (10x)	20 µl	180 µl	200 µl
Bead Wash Buffer (2x)	250 µl	250 µl	500 µl

2. **Vortex** each component and spin.
3. Transfer 100 µl of **1x Wash Buffer 1** into a fresh 0.2 ml PCR tube and pre-heat at **65°C**
4. Split the **1x Stringent Wash Buffer** into **two 200 µl aliquots, in 0.2 ml PCR tubes**, pre-heat at **65°C**.
5. Store the 200 µl of **1x Wash Buffer 1** and the **remaining 1x wash buffers at RT**.

Preparation of Dynabeads® M-270 Streptavidin

6. Mix **RT Dynabeads® M-270 Streptavidin** by vortexing for **15 seconds**.
7. Aliquot **50 µl** of Dynabeads® M-270 Streptavidin per capture reaction into a 1.5 ml tube.
8. Place the 1.5 ml tube on a magnet and incubate **20-30 seconds**.
9. Remove and discard the **supernatant**.
10. Add 200 µl of **1x Bead Wash Buffer** per capture reaction and vortex for **10 seconds**.
11. Repeat steps **8-10** once more for a total of **two washes**.
12. Place the 1.5 ml tube in a magnetic stand and incubate **20-30 seconds**.
13. Remove and discard the **supernatant**.
14. Add 100 µl of **1x Bead Wash Buffer** and vortex briefly.
15. Transfer 100 µl of resuspended beads into a new tube. **Proceed to the next step only when the hybridization incubation ends**.
16. Place the tube on a magnet and incubate **1-2 minutes** or until all beads have separated.
17. Remove and discard the **supernatant** and **proceed immediately to the next step**.




Protocol Guide v1.2.2

Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

Procedure

18. Set a thermocycler at **65°C to hold** with heated lid **at 70°C**.
19. Transfer the hybridization reaction mix to the tube containing the pelleted **Dynabeads® M-270 Streptavidin**.
20. Remove the 0.2 ml PCR tube from the magnet and mix the hybridization reaction mix with the Dynabeads® M-270 Streptavidin by pipetting **10 times**.
21. Transfer to the **65°C** thermocycler (heated lid at **70°C**), incubate for **45 minutes**.
22. **Every 12 minutes** during the 45-minute incubation at 65°C, vortex for **3 seconds**.
23. Remove the tube from the thermocycler, add **100 µl** of **pre-heated 1x Wash Buffer 1**.
24. Mix by pipetting **10 times** and place on a magnet for **10-15 seconds**.
25. Remove the **supernatant**.
26. Remove the tube from the magnet, add **200 µl** of pre-heated **1x Stringent Wash Buffer**
27. Mix by pipetting **10 times**.
28. Transfer to **65°C** thermocycler (heated lid at to **70°C**) and incubate **for 5 minutes**.
29. Remove from the thermocycler and place on a magnet stand for **10-15 seconds**.
30. Repeat steps 25-29 for a total of **two washes** with **heated 1x Stringent Wash Buffer**.
31. Remove the **supernatant**.
32. Remove from the magnet, add **200 µl** of **RT 1x Wash Buffer 1**.
33. **Vortex** for **2 minutes** and briefly spin.
34. Place the tube on a magnet for **20-30 seconds**.
35. Remove the **supernatant**.
36. Remove from the magnet, add **200 µl** of **RT 1x Wash Buffer 2**.
37. **Vortex** for **1 minute** and briefly spin.
38. Place the tube on a magnet for **20-30 seconds**.
39. Remove the **supernatant**.
40. Remove from the magnet, add **200 µl** of **RT 1x Wash Buffer 3**.
41. **Vortex** for **30 seconds** and briefly spin.
42. Place the tube on a magnet **for 1-2 minutes**.
43. Remove the **supernatant**.
44. Remove the tube from the magnet, resuspend in **24 µl** of **nuclease-free water** by pipetting **10-15 times**.



Protocol Guide v1.2.2

Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing

2.C Captured library amplification and clean-up

Before you start

- Thaw the PCR Master Mix – PostCap (2x) (green cap) and the Primer Mix – PostCap (10 µM) (**black cap**) **on ice**, vortex, spin down. Keep **on ice**.
- Equilibrate the **Target Pure™ NGS clean-up beads** to **RT** for **20-30 minutes**
- Prepare a solution of 80% Ethanol / 20% molecular biology grade water (400 µl per capture reaction)

Procedure

1. Set up the program.

Step	Temperature	Time	Cycles
1	98°C	Hold	1
2	98°C	2 min	1
3	98°C	20 sec	7-16*
4	60°C	30 sec	
5	72°C	30 sec	
6	72°C	1 min	1
7	4°C	Hold	1


** as a general guideline: 7 cycles for exome, 9 cycles for probe sets >5 Mb in size, 10 cycles for 3-5 Mb, 12 cycles for 0.2-3 Mb, 14 cycles for 0.2-0.04 Mb, 16 cycles for <0.04 Mb.*

Set the lid to heated at **105°C**, volume to **50 µl** and start.

2. Prepare the PCR reaction mix **on ice**. Mix by **pipetting 10 times** or **vortex** and **spin**.

Components	Volume for 1 reaction
PCR Master Mix – PostCap (2x)	25 µl
Primer Mix – PostCap (10 µM)	2.5 µl
Total	27.5 µl

3. Transfer **22.5 µl of resuspended Dynabeads® M-270 Streptavidin** with captured library DNA to the **27.5 µl** of PCR reaction master mix and mix by pipetting **10-15 times**.
4. Transfer to the pre-heated thermocycler (**98°C**) and skip to the next step.
5. Remove the samples from the cycling block and **proceed immediately to library clean-up using Target Pure™ NGS clean-up beads**.



Cell3™ Target: DNA Target Enrichment for Next Generation Sequencing


Clean-up of amplified captured library

6. Add **75 µl** of vortexed **RT Target Pure™ NGS clean-up beads** to a fresh tube.
7. Transfer the whole **50 µl** of PCR product to the 75 µl of Target Pure™ NGS clean-up beads and mix by **pipetting 15-20 times**, taking care to avoid the formation of bubbles.
8. **Incubate** the mixture for **5 minutes** at **RT**.
9. Place the tube on the magnet for **5 minutes** at **RT**.
10. Keeping the tube on the magnet, remove and **discard** the **supernatant**.
11. Add 200 µl of 80% ethanol to the tube at **RT** for **30 seconds**.
12. Repeat steps 10-11 for a total of **two** 80% ethanol washes.
13. Keeping the tube on the magnet, remove and **discard** the **supernatant**.
14. Use a **10 µl pipette** to remove any residual liquid.
15. Keeping the tube on the magnet, incubate at **RT** with **open lids** for **5 minutes** or until the beads are dry. **Avoid over-drying of beads.**
16. Remove the tube from the magnet and resuspend the dried beads in **32.5 µl** of **Buffer EB** or equivalent pipetting **10-15 times**.
17. **Incubate** the tube for **2 minutes** at **RT**.
18. Place the tube on the magnet for **2 minutes** at **RT**.
19. **Recover 30 µl** of **supernatant** and transfer to a fresh 1.5 ml **low-bind tube**.


Stop point: the captured DNA library can be stored at **-20°C**, if not proceeding immediately to the library quality check step.

2.D Captured library QC

- DNA **concentration** (ng/µl) and total yield (ng) (Qubit HS)
- DNA **average** fragment size and absence of additional lower or higher molecular weight peaks (TapeStation HS)



Experienced User



www.nonacus.com

Nonacus Limited

Birmingham Research Park, Vincent Drive, Birmingham, B15 2SQ

info@nonacus.com

© 2020 Nonacus Limited. For research use only. Not for use in diagnostic procedures.

NONACUS
Advancing Non-Invasive Healthcare