

Cell3™ Preserver – whole blood stabilization tube

For collection and stabilization of whole blood for circulating cell-free DNA (cfDNA) applications

Highlights

Stabilizes whole blood for up to 15 days at room temperature

Cell stability for up to 15 days post blood draw make Cell3 Preserver tubes ideal for batching samples prior to countrywide or international shipping.

Plastic evacuated tube with standard tube format

Avoids potential glass breakage and ensures ease of use by phlebotomists using standard draw equipment and disposables.

Formatted for clinical high throughput and automated workflows

Barcoded labels allow seamless workflow integration with automation, without the need to add barcodes manually.

Validated for non-invasive applications: NIPT and liquid biopsy testing

Quality and quantity of cfDNA extracted from blood and collected in the Cell3 Preserver tube was assessed for NIPT and liquid biopsy applications alongside the market leading cell stabilization tube to ensure accuracy of genetic testing.

Data generated across multiple platforms

Data generated using qPCR, DNA fragment analyzer and NGS (whole genome sequencing and targeted resequencing approaches) showed equivalent or better results/metrics for cfDNA analysis compared to the market leading company.

Introduction

Modern genetic techniques such as next generation sequencing and droplet digital PCR are extremely sensitive and can provide a rapid and accurate assessment of genetic variants from a blood sample. Recently, circulating cell-free (cf) DNA in biological fluids has been of growing interest in oncology, prenatal and transplant studies. However, blood collection in routine EDTA tubes can have a negative impact on cfDNA based tests if plasma is not isolated within 24 hours from blood draw, due to the degradation of nucleated (white blood) cells and shedding of genomic DNA into the plasma. Therefore, any delays in transport of samples from collection site to analysis location can negatively impact on results.

Tube format

Cell3 Preserver tubes consist of sterile, evacuated, polyethylene terephthalate tubes that are designed for direct-draw blood collection. They contain an additive at the correct volume to simultaneously stabilize blood cells and impede coagulation at the time of collection. The stabilizer acts by preserving the nucleated cells (white blood cells) in an intact state until processing and analysis can be performed. Phlebotomists can use these tubes with standard draw techniques and equipment. Tubes are labelled with barcodes which can be linked to patient details for LIMS compatibility and sample tracking to assist with high throughput and automated procedures.

Data quality

For quality assurance purposes, Cell3 Preserver tubes were compared to Streck tubes in key ccfDNA testing applications, such as NIPT and liquid biopsy.

NIPT: Serenity NIPT test (Illumina technology, based on the Verifi test) kindly undertaken at CooperGenomics. Blood samples from 5 patients were collected into both Cell3 Preserver and Streck tubes for testing. Plasma was isolated at 5 days after blood draw for Streck tubes and 5 and 11 days for Cell3 Preserver tubes. cfDNA was extracted and analyzed with the Serenity test for NIPT of fetal aneuploidies. Correct outcomes were observed for all tubes and time points (ie, 2 male and 2 female euploid, and one monosomy X). Testing metrics showed comparable data across all time points (Figures 1 and 2).




Figure 1: DNA concentration of PCR-free libraries prepared from ccfDNA for whole genome sequencing on NextSeq500 (Illumina).




Figure 2: Fetal fraction observed

Liquid biopsy: circulating tumor DNA (ctDNA) is showing great promise as a biomarker analyte for cancer diagnosis, stratification of patient treatment and disease surveillance. To validate Cell3 Preserver technology in this field, a blood sample from a patient with metastatic breast cancer was collected in Cell3 Preserver and Streck tubes post chemotherapy.

ctDNA was extracted and tested alongside gDNA obtained from leukocytes and a FFPE biopsy sample of the primary tumor. Targeted resequencing across a panel of genes commonly associated with cancer was performed at high depth using Illumina sequencing. Bioinformatics analysis was conducted to annotate variants and identify causative actionable mutations. Variant frequencies were calculated and compared between samples collected in Cell3 Preserver and Streck tubes (Figure 3).




Figure 3: Two identified pathogenic variants were found present in the FFPE and ctDNA samples, but not in the germline gDNA.

Summary

The Cell3 Preserver is a plastic, whole blood preservation tube validated for NIPT and liquid biopsy applications. Showing superior blood cell stability compared with the leading market company at over 13 days after blood draw, the Cell3 Preserver tube has been designed for sample collection and shipment from around the world to centralized lab facilities. A sample barcode is included to allow LIMS compatibility, sample tracking and to assist high throughput clinical and translational labs.

Learn more

To learn more about the Cell3 Preserver tubes and to download the protocols, please visit: www.nonacus.com

Ordering information

Product

Cell3 Preserver 9 ml, 50 Tube Pack

Catalog No.

PRE_C3P_WBS_50

Nonacus Limited
Quinton Business Park
11 Ridgeway
Birmingham
B32 1AF

info@nonacus.com

© 2023 Nonacus Limited. For research use only. Not for use in diagnostic procedures.